

Bard Woven

diagnostic catheter

BARD

Diagnostic performance without compromises

After more than 30 years, the Bard Woven diagnostic catheter is still setting the standard for reliability and performance. A workhorse in the truest sense, the Woven catheter has been used in more procedures than all plastic diagnostic catheters combined. Reliability, durable construction and clinical versatility have contributed to the Woven catheter's popularity, making it among the most trusted and successful catheters ever introduced.

Platinum electrodes connected to pure silver core conducting wires—used only by Bard—produce the highest resolution electrograms.

The Woven catheter's bioflexibility allows the shaft to soften rapidly at body temperature to a predetermined biosoftness.


- Biosoftening characteristic allows the catheter tip to conform to the endocardial wall for excellent tissue contact.
- Softening response prevents endocardial bounce, providing greater stability and better signal quality.
- Mapping accuracy is enhanced through the combination of bioflexibility, stable positioning, and high-resolution electrograms.



To produce woven bioflexibility, the woven core is impregnated with a proprietary polymer in an environment of strictly-controlled temperature and pressure.



Interweaving 16 strands of pure silver core conducting wire provides excellent conduction with low impedance.



More than 100 configurations are available to meet virtually any diagnostic need.

Ordering Information

Bard Woven Diagnostic Catheters

Item No.	French	Curve	Electrode Spacing	Poles	Usable Length
191001	6F	Cournand®	1cm	Bipolar	125cm
004019P	4F	Damato	1cm	Tripolar	105cm
004020P	5F	Damato	1cm	Tripolar	105cm
2001123	6F	Cournand	1cm	Tripolar	110cm
002854P	6F	Damato	1cm	Tripolar	105cm
200070	6F	Damato	1cm	Tripolar	110cm
200573	6F	Josephson	1cm	Tripolar	110cm
002853P	7F	Josephson	1cm	Tripolar	105cm
099021P	7F	Cournand Einstein	1cm	Tripolar	105cm
200071	7F	Damato	1cm	Tripolar	110cm
200150	4F	Cournand	2.5,2mm	Quadpolar	80cm
200151	4F	Josephson	2.5,2mm	Quadpolar	80cm
200152	4F	Cournand	1mm	Quadpolar	80cm
200153	4F	Josephson	1mm	Quadpolar	80cm
200055P	5F	Cournand	1cm	Quadpolar	110cm
200056	5F	Cournand	1cm	Quadpolar	110cm
200056E	5F	Cournand	1cm	Quadpolar	110cm
200057P	5F	Josephson	1cm	Quadpolar	110cm
200058	5F	Josephson	1cm	Quadpolar	110cm
200058E	5F	Josephson	1cm	Quadpolar	110cm
200059P	5F	Josephson	5mm	Quadpolar	110cm
200060	5F	Josephson	5mm	Quadpolar	110cm
200060E	5F	Josephson	5mm	Quadpolar	110cm
200061P	5F	Straight	1cm	Quadpolar	110cm
200062E	5F	Straight	1cm	Quadpolar	110cm
200063P	5F	Damato	1cm	Quadpolar	110cm
200064P	5F	Josephson	2.5,2mm	Quadpolar	110cm
200577	5F	Josephson	5mm	Quadpolar	120cm
200577E	5F	Josephson	5mm	Quadpolar	120cm
200580	5F	Cournand	5mm	Quadpolar	110cm
200580E	5F	Cournand	5mm	Quadpolar	110cm
200581	5F	Cournand	2.5,2mm	Quadpolar	110cm
200581E	5F	Cournand	2.5,2mm	Quadpolar	110cm
200587	5F	Cournand	5mm	Quadpolar	110cm
200587E	5F	Cournand	5mm	Quadpolar	110cm
200588	5F	Cournand	2.5,2mm	Quadpolar	110cm
200588E	5F	Cournand	2.5,2mm	Quadpolar	110cm
201400	5F	Damato	1cm	Quadpolar	110cm
008567	6F	Josephson	5mm	Quadpolar	120cm
002943P	6F	Straight	1cm	Quadpolar	110cm
006245P	6F	Cournand	1cm	Quadpolar	110cm
006692P	6F	Damato	1cm	Quadpolar	110cm
200066	6F	Cournand	5mm	Quadpolar	125cm
200067	6F	Cournand	1cm	Quadpolar	110cm
200068	6F	Straight	1cm	Quadpolar	110cm
200069	6F	Damato	1cm	Quadpolar	110cm
200201P	6F	Josephson	1cm	Quadpolar	125cm
200202	6F	Josephson	1cm	Quadpolar	125cm
200472	6F	Cournand	2mm	Quadpolar	115cm
200473	6F	Josephson	2mm	Quadpolar	115cm
200570P	6F	Cournand	2.5,2mm	Quadpolar	125cm

Item No.	French	Curve	Electrode Spacing	Poles	Usable Length
200571P	6F	Cournand	5mm	Quadpolar	125cm
200572P	6F	Damato	5mm	Quadpolar	125cm
200574	6F	Cournand	5mm	Quadpolar	125cm
200578	6F	Damato	5mm	Quadpolar	110cm
200582	6F	Cournand	2.5,2mm	Quadpolar	110cm
200583	6F	Cournand	2.5,2mm	Quadpolar	125cm
200584	6F	Josephson	2.5,2mm	Quadpolar	125cm
200585	6F	Damato	2.5,2mm	Quadpolar	125cm
200586	6F	Damato	2.5,2mm	Quadpolar	125cm
200589	6F	Cournand	2.5,2mm	Quadpolar	110cm
200590	6F	Cournand	2.5,2mm	Quadpolar	125cm
200591	6F	Josephson	2.5,2mm	Quadpolar	125cm
200592	6F	Damato	2.5,2mm	Quadpolar	125cm
200593	6F	Damato	2.5,2mm	Quadpolar	125cm
200624	6F	Josephson	5mm	Quadpolar	120cm
200719	6F	Damato	1cm	Quadpolar	110cm
200728	6F	Straight	5mm	Quadpolar	120cm
200792	6F	Damato	5mm	Quadpolar	110cm
002944P	7F	Straight	1cm	Quadpolar	110cm
200203	7F	Damato	1cm	Quadpolar	125cm
200793	7F	Damato	5mm	Quadpolar	110cm
200594	6F	Josephson	5.5,5,170 from tip	Pentapolar	120cm
201401	6F	Josephson	5.5,5,282 from tip	Pentapolar	120cm

Bard Woven Multipolar Mapping Catheters

Item No.	French	Curve	Electrode Spacing	Poles	Usable Length
200204	5F	Cournand	5mm	Hexapolar	110cm
200204E	5F	Cournand	5mm	Hexapolar	110cm
200205	5F	Cournand	1cm	Hexapolar	110cm
200205E	5F	Cournand	1cm	Hexapolar	110cm
200579P	5F	Cournand	2.5,2mm	Hexapolar	100cm
200085P	6F	Straight	10,10,10,70,10mm	Hexapolar	110cm
200206	6F	Cournand	1cm	Hexapolar	125cm
200207	6F	Josephson	2.5,2mm	Hexapolar	125cm
200208	6F	Cournand	5mm	Hexapolar	125cm
200575	6F	Straight	2,11,2,8,3,2mm	Hexapolar	110cm
004178P	7F	Straight	10,80,10,20,10mm	Hexapolar	125cm
200209	6F	Cournand	1mm	Octapolar	125cm
200210	6F	Cournand	5mm	Octapolar	125cm
200354	6F	Cournand	1,10,1mm	Octapolar	125cm
200567	5F	Small 60	2.5,2mm	Decapolar	65cm+
200568	5F	Small 60	2.8,2mm	Decapolar	65cm+
200211	6F	Cournand	2.5,2mm	Decapolar	125cm
200212	6F	Cournand	2mm	Decapolar	125cm
200355	6F	Cournand	5mm	Decapolar	125cm
200576	6F	Damato	2.5,2mm	Decapolar	110cm
200595	6F	Small 60	2.5,2mm	Decapolar	110cm
200612P	6F	Josephson	5.5,5,6,10,10,10,10mm	Decapolar	125cm

Bard Woven Open Lumen Catheters

Item No.	French	Curve	Electrode Spacing	Poles	Usable Length
200076P	6F	Cournand	1cm	Quadpolar	110cm
200078P	6F	Damato	1cm	Quadpolar	110cm
200077P	6F	Josephson	5mm	Quadpolar	110cm

Bard Woven Cable Information

Item No.	Use	Length	Description
008568P	Single Use	125cm	Protected 4-Pin cable for use with Woven
200088P	Single Use	125cm	Protected 4-Pin Easy-Mate™ for use with Woven
200089P	Single Use	125cm	Protected 6-Pin Easy-Mate for use with Woven
560004P	Single Use	120cm	Protected 10-Pin SureLink™ for use with Woven+
560001P	Single Use	210cm	Protected 10-Pin SureLink for use with Woven+
560004A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	120cm	Protected 10-Pin SureLink for use with Woven+
560001A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	210cm	Protected 10-Pin SureLink for use with Woven+

*Not available in all countries. Please check with your local Bard representative for availability.

Bard Electrophysiology Division
USA
Tel: 978-441-6202
800-824-8724 (toll free w/in the US)
Fax: 978-323-2222
877-281-1492 (toll free w/in the US)

Bard Asia Pacific Division
USA
Tel: 908-277-8000
Fax: 908-277-8090

Bard Australia, Pty Ltd
AUSTRALIA
Tel: 61-2-8875-4000
Fax: 61-2-8875-4050

Bard Benelux, N.V.
BELGIUM
Tel: 32-0-1428-5900
Fax: 32-0-1428-6955

Bard Canada, Inc
CANADA
Tel: 905-275-8000
Fax: 905-275-5370

Bard International, Inc
CHINA
Tel: 86-10-6501-4225
Fax: 86-10-6501-2573

Laboratories Bard SA
FRANCE
Tel: 33-1-39-30-5858
Fax: 33-1-39-30-5859

C. R. Bard Germany, GmbH
GERMANY
Tel: 49-721-9445-265
Fax: 49-721-9445-111

Bard International
INDIA
Tel: 91-22-612-3783
Fax: 91-22-612-3821

Bard S.p.A.
ITALY
Tel: 39-06-52-493-1
Fax: 39-06-52-95067

Medicon, Inc
JAPAN
Tel: 6-6203-6546
Fax: 6-6203-1516

Bard International, Inc.
KOREA
Tel: 82-2-719-3080
Fax: 82-2-719-3081

Bard Latin America
USA
Tel: 908-277-8000
Fax: 908-277-8090

C. R. Bard Portugal, LDA
PORTUGAL
Tel: 351 21-319-0330
Fax: 351 21-319-0339

Bard Singapore Pte. Ltd.
SINGAPORE
Tel: 65-250-7311
Fax: 65-353-7311

Bard De Espana, SA
SPAIN
Tel: 34-93-253-7800
Fax: 34-93-205-8229

Bard Medical Systems Norden
SWEDEN
Tel: 46-42-386-000
Fax: 46-42-386-010

Bard Medica, SA
SWITZERLAND
Tel: 41-17-225-360
Fax: 41-17-225-370

Pacific Health Care Co., Ltd
TAIWAN
Tel: 886-2274-07-891
Fax: 886-2277-75-576

Bard Limited
UNITED KINGDOM
Tel: 44-1-293-527-888
Fax: 44-1-293-552-428

Bard Sdn Bhd
WEST MALAYSIA
Tel: 60-4-489-1300
Fax: 60-4-489-1737



Bard Electrophysiology
55 Technology Drive • Lowell, MA 01851 USA
Tel: (978) 441-6202 • Fax: (978) 323-2222
www.BARDEP.com



Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and use.

Bard and the stylized heart and Cournand are registered trademarks of C.R. Bard, Inc. or an affiliate. Josephson, Easy-Mate and SureLink are trademarks of C.R. Bard, Inc. or an affiliate. All trademarks or registered trademarks are the property of their respective owners.
© 2002 C.R. Bard, Inc. All rights reserved. LT0420003/Rev. 01/12-02