

Bard EP•XT™ steerable catheter

BARD

Incredible accuracy and control to speed up and simplify intracardiac mapping

The reliable, easy-to-steer Bard EP•XT

catheter offers an exceptional blend of shaft stability

and curve control for faster, more accurate mapping.

And Bard's exclusive point-of-tip actuation control

supplies exceptional levels of tip response and control.

This extraordinary family of small, medium and long-

reach curves provides clinical versatility and reliable

performance in a wide variety of human anatomies.

Gentle curve shape enables you to meet individual anatomical challenges.

The spring-coil core and wire-braid outer jacket combine to provide superior column strength.

Soft, durometer polyurethane eliminates distal tip trauma.

The EP•XT catheter is available in 6F and 7F shafts, as well as your choice of 2mm and 4mm distal tip electrodes.

Patented curving mechanism allows smooth, precise incremental movement.



Curve shape control technology prevents unwanted tip movement.



Convenient Easy-Mate™ connectors make cabling up fast and easy.

Unique, user-friendly thumbwheel handle provides exceptional control.



Ordering Information

Bard EP•XT™ Steerable Catheters

Item No.	French	Electrode Curve	Spacing	Poles	Tip Electrode	Proximal Electrode	Usable Length
200769	6F	Large	5mm	Quadpolar	2mm	1mm	110cm
200770	6F	Large	2,5,2mm	Quadpolar	2mm	1mm	110cm
201116	6F	Large	1mm	Quadpolar	2mm	1mm	110cm
200771	6F	Large	1cm	Hexapolar	2mm	1mm	110cm
200772	6F	Large	2,5,2mm	Hexapolar	2mm	1mm	110cm
201009	6F	Large	5mm	Hexapolar	2mm	1mm	110cm
200794	6F	Large	2,5,2mm	Octapolar	2mm	1mm	110cm
200795	6F	Large	2mm	Octapolar	2mm	1mm	110cm
200796	6F	Large	5mm	Octapolar	2mm	1mm	110cm
200797	6F	Large	2,10,2mm	Octapolar	2mm	1mm	110cm
200798	6F	Large	2,1,1,1,1,1mm	Octapolar	2mm	1mm	110cm
201007	6F	Large	2,5,2mm	Decapolar	2mm	1mm	110cm

Bard EP•XT Cable Information

Item No.	Use	Length	Description
200088P	Single Use	125cm	Protected 4-Pin Easy-Mate™ for use with EP•XT Quadpolar
200089P	Single Use	125cm	Protected 6-Pin Easy-Mate for use with EP•XT Hexapolar
200774P	Single Use	125cm	Protected 8-Pin Easy-Mate for use with EP•XT Octapolar
560004P	Single Use	120cm	Protected 10-Pin SureLink™ for use with EP•XT Decapolar
560001P	Single Use	210cm	Protected 10-Pin SureLink for use with EP•XT Decapolar
560004A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	120cm	Protected 10-Pin SureLink for use with EP•XT Decapolar
560001A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	210cm	Protected 10-Pin SureLink for use with EP•XT Decapolar

*Not available in all countries. Please check with your local Bard representative for availability.

Bard Electrophysiology Division
USA
Tel: 978-441-6202
800-824-8724 (toll free w/in the US)
Fax: 978-323-2222
877-281-1492 (toll free w/in the US)

Bard Asia Pacific Division
USA
Tel: 908-277-8000
Fax: 908-277-8090

Bard Australia, Pty Ltd
AUSTRALIA
Tel: 61-2-8875-4000
Fax: 61-2-8875-4050

Bard Benelux, N.V.
BELGIUM
Tel: 32-0-1428-5900
Fax: 32-0-1428-6955

Bard Canada, Inc
CANADA
Tel: 905-275-8000
Fax: 905-275-5370

Bard International, Inc
CHINA
Tel: 86-10-6501-4225
Fax: 86-10-6501-2573

Laboratories Bard SA
FRANCE
Tel: 33-1-39-30-5858
Fax: 33-1-39-30-5859

C. R. Bard Germany, GmbH
GERMANY
Tel: 49-721-9445-265
Fax: 49-721-9445-111

Bard International
INDIA
Tel: 91-22-612-3783
Fax: 91-22-612-3821

Bard S.p.A.
ITALY
Tel: 39-06-52-493-1
Fax: 39-06-52-95067

Medicon, Inc
JAPAN
Tel: 6-6203-6546
Fax: 6-6203-1516

Bard International, Inc.
KOREA
Tel: 82-2-719-3080
Fax: 82-2-719-3081

Bard Latin America
USA
Tel: 908-277-8000
Fax: 908-277-8090

C. R. Bard Portugal, LDA
PORTUGAL
Tel: 351 21-319-0330
Fax: 351 21-319-0339

Bard Singapore Pte. Ltd.
SINGAPORE
Tel: 65-250-7311
Fax: 65-353-7311

Bard De Espana, SA
SPAIN
Tel: 34-93-253-7800
Fax: 34-93-205-8229

Bard Medical Systems Norden
SWEDEN
Tel: 46-42-386-000
Fax: 46-42-386-010

Bard Medica, SA
SWITZERLAND
Tel: 41-17-225-360
Fax: 41-17-225-370

Pacific Health Care Co., Ltd
TAIWAN
Tel: 886-2274-07-891
Fax: 886-2277-75-576

Bard Limited
UNITED KINGDOM
Tel: 44-1-293-527-888
Fax: 44-1-293-552-428

Bard Sdn Bhd
WEST MALAYSIA
Tel: 60-4-489-1300
Fax: 60-4-489-1737



Bard Electrophysiology
55 Technology Drive • Lowell, MA 01851 USA
Tel: (978) 441-6202 • Fax: (978) 323-2222
www.BARDEP.com



Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and use.

Bard and the stylized heart are registered trademarks of C.R. Bard, Inc. or an affiliate. EP•XT, Easy-Mate and SureLink are trademarks of C.R. Bard, Inc. or an affiliate. All trademarks or registered trademarks are the property of their respective owners.
© 2002 C.R. Bard, Inc. All rights reserved. LT0420008/Rev. 01/12-02