

# Bard EP•XT™

steerable catheter

BARD

## Incredible accuracy and control to speed up and simplify intracardiac mapping

The reliable, easy-to-steer Bard EP•XT

catheter offers an exceptional blend of shaft stability

and curve control for faster, more accurate mapping.

And Bard's exclusive point-of-tip actuation control

supplies exceptional levels of tip response and control.

This extraordinary family of small, medium and long-

reach curves provides clinical versatility and reliable

performance in a wide variety of human anatomies.


Gentle curve shape enables you to meet individual anatomical challenges.

The spring-coil core and wire-braid outer jacket combine to provide superior column strength.


Soft, durometer polyurethane eliminates distal tip trauma.

The EP•XT catheter is available in 6F and 7F shafts, as well as your choice of 2mm and 4mm distal tip electrodes.

Patented curving mechanism allows smooth, precise incremental movement.




Curve shape control technology prevents unwanted tip movement.



Convenient Easy-Mate™ connectors make cabling up fast and easy.



Unique, user-friendly thumbwheel handle provides exceptional control.



# Ordering Information

## Bard EP•XT™ Steerable Catheters

Item No.	French	Electrode Curve	Spacing	Poles	Tip Electrode	Proximal Electrode	Usable Length
200769	6F	Large	5mm	Quadpolar	2mm	1mm	110cm
200770	6F	Large	2,5,2mm	Quadpolar	2mm	1mm	110cm
201116	6F	Large	1mm	Quadpolar	2mm	1mm	110cm
200771	6F	Large	1cm	Hexapolar	2mm	1mm	110cm
200772	6F	Large	2,5,2mm	Hexapolar	2mm	1mm	110cm
201009	6F	Large	5mm	Hexapolar	2mm	1mm	110cm
200794	6F	Large	2,5,2mm	Octapolar	2mm	1mm	110cm
200795	6F	Large	2mm	Octapolar	2mm	1mm	110cm
200796	6F	Large	5mm	Octapolar	2mm	1mm	110cm
200797	6F	Large	2,10,2mm	Octapolar	2mm	1mm	110cm
200798	6F	Large	2,1,1,1,1,1mm	Octapolar	2mm	1mm	110cm
201007	6F	Large	2,5,2mm	Decapolar	2mm	1mm	110cm

## Bard EP•XT Cable Information

Item No.	Use	Length	Description
200088P	Single Use	125cm	Protected 4-Pin Easy-Mate™ for use with EP•XT Quadpolar
200089P	Single Use	125cm	Protected 6-Pin Easy-Mate™ for use with EP•XT Hexapolar
200774P	Single Use	125cm	Protected 8-Pin Easy-Mate™ for use with EP•XT Octapolar
560004P	Single Use	120cm	Protected 10-Pin SureLink™ for use with EP•XT Decapolar
560001P	Single Use	210cm	Protected 10-Pin SureLink™ for use with EP•XT Decapolar
560004A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	120cm	Protected 10-Pin SureLink™ for use with EP•XT Decapolar
560001A*	Autoclave or Sterrad 10x Reuse or ETO 3x Reuse	210cm	Protected 10-Pin SureLink™ for use with EP•XT Decapolar

\*Not available in all countries. Please check with your local Bard representative for availability.

**Bard Electrophysiology Division**  
USA  
Tel: 978-441-6202  
800-824-8724 (toll free w/in the US)  
Fax: 978-323-2222  
877-281-1492 (toll free w/in the US)

**Bard Asia Pacific Division**  
USA  
Tel: 908-277-8000  
Fax: 908-277-8090

**Bard Australia, Pty Ltd**  
AUSTRALIA  
Tel: 61-2-8875-4000  
Fax: 61-2-8875-4050

**Bard Benelux, N.V.**  
BELGIUM  
Tel: 32-0-1428-5900  
Fax: 32-0-1428-6955

**Bard Canada, Inc**  
CANADA  
Tel: 905-275-8000  
Fax: 905-275-5370

**Bard International, Inc**  
CHINA  
Tel: 86-10-6501-4225  
Fax: 86-10-6501-2573

**Laboratories Bard SA**  
FRANCE  
Tel: 33-1-39-30-5858  
Fax: 33-1-39-30-5859

**C. R. Bard Germany, GmbH**  
GERMANY  
Tel: 49-721-9445-265  
Fax: 49-721-9445-111

**Bard International**  
INDIA  
Tel: 91-22-612-3783  
Fax: 91-22-612-3821

**Bard S.p.A.**  
ITALY  
Tel: 39-06-52-493-1  
Fax: 39-06-52-95067

**Medicon, Inc**  
JAPAN  
Tel: 6-6203-6546  
Fax: 6-6203-1516

**Bard International, Inc.**  
KOREA  
Tel: 82-2-719-3080  
Fax: 82-2-719-3081

**Bard Latin America**  
USA  
Tel: 908-277-8000  
Fax: 908-277-8090

**C. R. Bard Portugal, LDA**  
PORTUGAL  
Tel: 351 21-319-0330  
Fax: 351 21-319-0339

**Bard Singapore Pte. Ltd.**  
SINGAPORE  
Tel: 65-250-7311  
Fax: 65-353-7311

**Bard De Espana, SA**  
SPAIN  
Tel: 34-93-253-7800  
Fax: 34-93-205-8229

**Bard Medical Systems Norden**  
SWEDEN  
Tel: 46-42-386-000  
Fax: 46-42-386-010

**Bard Medica, SA**  
SWITZERLAND  
Tel: 41-17-225-360  
Fax: 41-17-225-370

**Pacific Health Care Co., Ltd**  
TAIWAN  
Tel: 886-2274-07-891  
Fax: 886-2277-75-576

**Bard Limited**  
UNITED KINGDOM  
Tel: 44-1-293-527-888  
Fax: 44-1-293-552-428

**Bard Sdn Bhd**  
WEST MALAYSIA  
Tel: 60-4-489-1300  
Fax: 60-4-489-1737



Bard Electrophysiology  
55 Technology Drive • Lowell, MA 01851 USA  
Tel: (978) 441-6202 • Fax: (978) 323-2222  
[www.BARDEP.com](http://www.BARDEP.com)



Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and use.

Bard and the stylized heart are registered trademarks of C.R. Bard, Inc. or an affiliate. EP•XT, Easy-Mate and SureLink are trademarks of C.R. Bard, Inc. or an affiliate. All trademarks or registered trademarks are the property of their respective owners.  
© 2010 C.R. Bard, Inc. All rights reserved. LT04Z0008/Rev. 02/12-2010