



Performance in Motion.

tango™


STABILENE™ CATHETER

# tango

STABILENE CATHETER


## Performance in

Other catheters bounce...



...the Bard Tango catheter hugs the endocardium

Stabilene polymer construction assures consistent, long-term wall contact with less repositioning.



Single braid

Bard Tango catheter's double-stranded braid

Bard Tango Stabilene catheter's double-stranded braid offers more torque control than other catheters.


The Tango™...a dance of touch. It's the ultimate combination of precision and flexibility, one that relies on seamless contact and synergy between partners. The same principle has led to the development of Bard's latest line of catheters, The Tango Stabilene™ Catheter.

### A Partnership in Harmony

The Tango line of catheters is designed to synchronize with the moving endocardium and to maintain contact throughout even the longest procedure. The dampening characteristic of the Stabilene polymer significantly reduces catheter bounce resulting in improved stability. It's no secret...better contact yields optimal signal quality and procedural efficiency.

### Strength and Elegance

The ideal catheter has always been a blend of both stability and subtlety and Tango was designed with that premise in mind. Tango combines an exclusive double-strand wire braid for enhanced column support, with a soft tip to strike the perfect balance between strength and elegance.



The Tango catheter  
combines a firm shaft  
with a soft, flexible tip.

Proprietary Stabilene polymer  
adds stability and high  
dampening characteristics.

# Motion.

## Precision and Control

Like a dance partner that can intuitively follow the lead, Tango responds to the physicians' slightest movement. Its unique double-end braid optimizes torque response for ultimate control and precision positioning.

## A Passion for Performance

Like those that embrace the Tango, Bard Electrophysiology has always held an unwavering commitment to performance. The Tango catheter embodies over 25 years of development experience. Combining proprietary materials and advanced construction techniques, Tango maximizes performance while meeting your need for a low profile catheter.

The Tango Stabilene Catheter Features:

- Proprietary Stabilene high dampening polymer for continual endocardial contact and improved stability
- Enhanced column strength with soft-tip construction
- Double-stranded braid for optimal torque transmission



# Ordering Information

## Tango Stabilene Catheter

Item #	French	Curve	Poles	Spacing	Length
400540	5F	Cournand	Bipolar	1cm	110cm
400539	5F	Josephson	Bipolar	1cm	110cm
400551	5F	Cournand	Quadpolar	1cm	110cm
400553	5F	Cournand	Quadpolar	2,5mm	110cm
400545	5F	Cournand	Quadpolar	5mm	110cm
400542	5F	Cournand	Quadpolar	2mm	110cm
400546	5F	Damato	Quadpolar	5mm	110cm
405123	5F	Hisser	Quadpolar	2,5mm	110cm
405122	5F	Hisser	Quadpolar	2mm	110cm
405124	5F	Hisser	Quadpolar	5mm	110cm
400550	5F	Josephson	Quadpolar	1cm	110cm
400552	5F	Josephson	Quadpolar	2,5mm	110cm
400541	5F	Josephson	Quadpolar	2mm	110cm
400544	5F	Josephson	Quadpolar	5mm	110cm
400549	5F	Interian	Pentapolar	10,5,5,5mm	110cm
400557	5F	Cournand	Hexapolar	5mm	110cm
400571	5F	Cournand	Decapolar	2mm	110cm
400575	5F	Cournand	Decapolar	2,5mm	110cm
400577	5F	CS	Decapolar	2,5mm	110cm
400579	5F	CS	Decapolar	2,8mm	110cm
400572	5F	CS	Decapolar	2mm	110cm
400576	5F	Elmhurst	Decapolar	2,5mm	110cm
400570	5F	Josephson	Decapolar	2mm	110cm
400574	5F	Josephson	Decapolar	2,5mm	110cm

## SureLink™ Cables

Item #	Description	Use	Length
560002P	4-pin SureLink	ETO	120cm
560002BP	4-pin SureLink (Blue)	ETO	120cm
560002RP	4-pin SureLink (Red)	ETO	120cm
560002YP	4-pin SureLink (Yellow)	ETO	120cm
560003P	4-pin SureLink	ETO	210cm
560004P	10-pin SureLink	ETO	120cm
560001P	10-pin SureLink	ETO	210cm
560002A*	4-pin SureLink	Autoclave, Sterrad or ETO	120cm
560002BA*	4-pin SureLink (Blue)	Autoclave, Sterrad or ETO	120cm
560002RA*	4-pin SureLink (Red)	Autoclave, Sterrad or ETO	120cm
560002YA*	4-pin SureLink Yellow	Autoclave, Sterrad or ETO	120cm
560003A*	4-pin SureLink	Autoclave, Sterrad or ETO	210cm
560004A*	10-pin SureLink	Autoclave, Sterrad or ETO	120cm
560001A*	10-pin SureLink	Autoclave, Sterrad or ETO	210cm

\* Not available in all countries. Please check with your local Bard representative for availability.

For more information please call (978) 441-6202 or (800) 824-8724 or visit our website at [www.BardEP.com](http://www.BardEP.com)



Bard Electrophysiology Division  
 C. R. Bard, Inc.  
 55 Technology Drive • Lowell, MA 01851 USA  
 Tel: (978) 441-6202 • Fax: (978) 323-2222  
[www.BARDEP.com](http://www.BARDEP.com)



**Bard Electrophysiology Division**  
 USA  
 Tel 978-441-6202  
 800-824-8724 (toll free w/in the US)  
 Fax 978-323-2222  
 877-281-1492 (toll free w/in the US)

**Bard Asia Pacific Division**  
 USA  
 Tel 908-277-8000  
 Fax 908-277-8090  
**Bard Australia, Pty Ltd**  
 AUSTRALIA  
 Tel 61-2-8875-4000  
 Fax 61-2-8875-4050

**Bard Medica S.A.**  
 AUSTRIA  
 Tel 43-1-4949-130  
 Fax 43-1-4949-13030

**Bard Benelux, N.V.**  
 BELGIUM  
 Tel 32-0-1428-5900  
 Fax 32-0-1428-6955

**Bard Canada, Inc.**  
 CANADA  
 Tel 905-275-8000  
 Fax 905-275-5370

**Bard International, Inc.**  
 CHINA  
 Tel 86-10-6501-4225  
 Fax 86-10-6501-2573

**Laboratories Bard SA**  
 FRANCE  
 Tel 33-1-39-30-5858  
 Fax 33-1-39-30-5859

**C. R. Bard Germany, GmbH**  
 GERMANY  
 Tel 49-721-9445-265  
 Fax 49-721-9445-111

**Bard South Asia**  
 INDIA  
 Tel 91-22-612-3783  
 Fax 91-22-612-3821

**Bard S.p.A.**  
 ITALY  
 Tel 39-06-52-493-1  
 Fax 39-06-52-950-67

**Medicon, Inc.**  
 JAPAN  
 Tel 6-6203-6546  
 Fax 6-6203-1516

**Bard International, Inc.**  
 KOREA  
 Tel 82-2-719-3080  
 Fax 82-2-719-3081

**Bard Latin America**  
 USA  
 Tel 908-277-8000  
 Fax 908-277-8090

**C. R. Bard Portugal, LDA**  
 PORTUGAL  
 Tel 351-21-319-0330  
 Fax 351-21-319-0339

**Bard International, Inc.**  
 SINGAPORE  
 Tel 65-250-7311  
 Fax 65-353-7311

**Bard De Espana, SA**  
 SPAIN  
 Tel 34-93-253-7800  
 Fax 34-93-205-8229

**Bard Medical Systems Norden**  
 SWEDEN  
 Tel 46-42-386-000  
 Fax 46-42-386-010

**Bard Medica, SA**  
 SWITZERLAND  
 Tel 41-17-225-360  
 Fax 41-17-225-370

**Pacific Health Care Co., Ltd**  
 TAIWAN  
 Tel 886-2274-07-891  
 Fax 886-2277-75-576

**Bard Limited**  
 UNITED KINGDOM  
 Tel 44-1-293-527-888  
 Fax 44-1-293-552-428

**Bard Sdn Bhd**  
 WEST MALAYSIA  
 Tel 60-4-489-1300  
 Fax 60-4-489-1737

Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and use.

Bard and the stylized heart are registered trademarks of C.R. Bard, Inc. or an affiliate. Tango, Stabilene and SureLink are trademarks of C.R. Bard, Inc. or an affiliate. All trademarks or registered trademarks are the property of their respective owners.

© 2003 C.R. Bard, Inc. All rights reserved. LT04Z0044/Rev. 01/07-03